

Regular

Product Availability and Markings

| Material | Index | Power Range | Nominal Base Curve | Size | Add | Color | Fitting Height | Invisible Marking | |
|-----------------------------|-------|-----------------|---------------------|-------|--------------|----------------------|----------------|-------------------|----------|
| Hard Resin | 1.498 | -8.25 to +5.75 | 1.5, 3, 5, 6.5, 7.5 | 75/80 | 0.75 to 3.50 | Clear | 19mm | X | X 250 |
| Superlite™ 1.6 | 1.6 | -11.25 to +7.00 | 1.5, 3, 5, 6.5, 7.5 | 75/80 | 0.75 to 3.50 | Clear | 19mm | X H | X 250 |
| Superlite 1.6 Transitions® | 1.6 | -11.25 to +7.00 | 1.5, 3, 5, 6.5, 7.5 | 72/77 | 0.75 to 3.50 | Gray & Brown | 19mm | X H | X 250 |
| Superlite 1.67 | 1.67 | -12.50 to +7.50 | 1.5, 3, 5, 6.5, 7.5 | 75/80 | 0.75 to 3.50 | Clear | 19mm | ✧ H67 | ✧ 250 |
| Superlite 1.67 Transitions® | 1.67 | -12.50 to +7.50 | 1.5, 3, 5, 6.5, 7.5 | 72/77 | 0.75 to 3.50 | Gray & Brown | 19mm | ✧ H67 | ✧ 250 |
| Glass | 1.523 | -5.25 to +6.00 | 3,5,7 | 71/76 | 0.75 to 3.50 | Clear, PSSB and PSSG | 18mm | X | X 250 |
| Glass | 1.6 | -9.00 to +6.00 | 1.5, 3, 5, 7 | 70/75 | 0.75 to 3.50 | Clear, PGX, PBX | 19mm | X H | X 250 |
| Glass | 1.7 | -12.00 to -6.00 | 2 | 70/75 | 0.75 to 3.50 | Clear | 19mm | X H7 | X 250 |
| Glass | 1.8 | -20.0 to -4.25 | 1.5, 2 | 65/69 | 1.00 to 3.00 | Clear | 19mm | ◇ 18 | ◇ 25 |