

Short

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve	Size	Add	Color	Fitting Height	Invisible Marking	
Hard Resin	1.498	-8.50 to +5.75	1.5, 3, 5, 6.5, 7.5	75/80	0.75 to 3.50	Clear	16mm	○	○ 250
Transitions®	1.5	-8.50 to +6.00	1.5, 3, 5, 6.5, 7.5	72/77	0.75 to 3.50	Gray & Brown	16mm	○	○ 250
SuperLite™ 1.6	1.6	-11.25 to +7.00	1.5, 3, 5, 6.5, 7.5	75/80	0.75 to 3.50	Clear	16mm	○ H	○ 250
SuperLite 1.6 Transitions®	1.6	-11.25 to +7.00	1.5, 3, 5, 6.5, 7.5	72/77	0.75 to 3.50	Gray & Brown	16mm	○ H	○ 250
SuperLite 1.67	1.67	-12.50 to +6.00	1.5, 3, 4.5, 5.75, 7	75/80	1.00 to 3.00	Clear	16mm	○ H67	○ 250
SuperLite 1.67 Transitions®	1.67	-12.50 to +6.00	1.5, 3, 4.5, 5.75, 7	72/77	1.00 to 3.00	Gray & Brown	16mm	○ H67	○ 250
Glass	1.523	-5.25 to +6.00	3, 5, 7	71/76	0.75 to 3.50	Clear & PGX	16mm	○	○ 250
Glass	1.523	-5.25 to +6.00	3, 5, 7	71/76	0.75 to 3.50	T&D PGX	16mm	TD ○	○ 250
Glass	1.7	-12.00 to -6.00	2	70/75	0.75 to 3.50	Clear	16mm	⬡ 17	⬡ 250