

Product Availability and Markings

| Material | Index | Power Range | Base Curve | Size | Dynamic Reduction (D)* | Color | Fitting Height |
|------------|-------|----------------|------------|-------|------------------------|-------|----------------|
| Hard Resin | 1.498 | -6.00 to +6.00 | 2.5, 5, 7 | 75/80 | 0.75, 1.25, 1.75 | Clear | 16mm |
| Glass | 1.5 | -6.00 to +6.00 | 2.5, 5, 7 | 71/76 | 0.75, 1.25, 1.75 | Clear | 16mm |

* D.P. 2.25 (D) is not available in base 2.5