

Shamir Autograph Intelligence™

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Hight	Invisible Marking	
Hard Resin	1.498	-12.25 to +6.50	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	11 mm and up	15 AI	250
Blue Zero	1.498	-12.25 to +6.50	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	11 mm and up	BZ5 AI	250
Transitions® Signature®	1.5	-12.25 to +6.50	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray, Brown & Graphite Green	11 mm and up	TR AI	250
Transitions® XTRActive®	1.5	-12.25 to +6.50	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray	11 mm and up	X5 AI	250
Polarized	1.498	-12.25 to +6.50	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray & Brown	11 mm and up	PZ AI	250
Drivewear®	1.5	-12.25 to +6.50	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Drivewear®	11 mm and up	DW AI	250
DLC™ (Trivex)	1.53	-13.25 to +7.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	11 mm and up	TX AI	250
DLC™ Blue Zero	1.53	-13.25 to +7.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	11 mm and up	BZX AI	250
DLC™ (Trivex) Transitions® Signature®	1.53	-13.25 to +7.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray, Brown & Graphite Green	11 mm and up	TTX AI	250
DLC™ (Trivex) Transitions® XTRActive®	1.53	-13.25 to +7.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray	11 mm and up	XX AI	250
DLC™ (Trivex) Polarized	1.53	-13.25 to +7.00	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray & Brown	11 mm and up	PTX AI	250
1.56	1.56	-14.00 to +7.50	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	11 mm and up	56 AI	250

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.

Shamir Autograph Intelligence™

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Hight	Invisible Marking	
Polycarbonate	1.586	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	11 mm and up	PC S AI	S 250
Polycarbonate Blue Zero	1.586	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	11 mm and up	BZP S AI	S 250
Polycarbonate Transitions® Signature®	1.586	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray, Brown & Graphite Green	11 mm and up	PT S AI	S 250
Polycarbonate Transitions® XTRActive®	1.586	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray	11 mm and up	PX S AI	S 250
Polycarbonate Polarized	1.586	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray & Brown	11 mm and up	PPZ S AI	S 250
Polycarbonate Drivewear®	1.586	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Drivewear®	11 mm and up	PD S AI	S 250
SuperLite™ 1.6	1.6	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	11 mm and up	16 S AI	S 250
SuperLite™ 1.6 Blue Zero	1.6	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	11 mm and up	BZ6 S AI	S 250
SuperLite 1.6 Transitions® Signature®	1.6	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray, Brown & Graphite Green	11 mm and up	T6 S AI	S 250
SuperLite 1.6 Transitions® XTRActive®	1.6	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray	11 mm and up	X6 S AI	S 250
SuperLite 1.6 Polarized	1.6	-14.75 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray & Brown	11 mm and up	PZ6 S AI	S 250
SuperLite 1.67	1.67	-16.75 to +10.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	11 mm and up	67 S AI	S 250

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Hight	Invisible Marking
SuperLite™ 1.67 Blue Zero	1.67	-16.75 to +10.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	11 mm and up	BZ67 AI 250
SuperLite 1.67 Transitions® Signature®	1.67	-16.75 to +10.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray, Brown & Graphite Green	11 mm and up	T67 AI 250
SuperLite 1.67 Transitions® XTRActive®	1.67	-16.75 to +10.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray	11 mm and up	X67 AI 250
SuperLite 1.67 Polarized	1.67	-16.75 to +10.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray & Brown	11 mm and up	PZ67 AI 250
Superlite 1.74	1.74	-18.75 to +15.50	0.5, 1.5, 2.5, 3.5, 4.75, 6, 7.25, 8.25, 9, 9.75, 10.5, 11.25, 12	0.75 to 3.50	Clear	11 mm and up	74 AI 250
SuperLite 1.74 Transitions® Signature®	1.74	-18.75 to +15.50	0.5, 1.5, 2.5, 3.5, 4.75, 6, 7.25, 8.25, 9, 9.75, 10.5, 11.25, 12	0.75 to 3.50	Gray, Brown & Graphite Green	11 mm and up	T74 AI 250

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.